

## ProBasics Economy Two-Button Folding Steel Walker, (Adult)

HCPCS:

- ProBasics Economy Two-Button Folding Steel Walker provides maximum strength in a durable design, all backed by a limited lifetime warranty.
- Ample height adjustment options coupled with a two-button folding mechanism make it easy to use and accommodating to the most demanding patient needs.
- Available with or without 5" wheels and in junior or adult sizes.
- Supports patient weights up to 350 pounds.
- 1-inch steel frame construction provides superior strength
- Easy-to-use push button folding mechanism provides audible click when locked
- Each side operates independently to allow easy movement through narrow spaces and greater stability while standing
- Front cross brace near top of walker allows for full stride
- Side bar adjusts with height to provide additional support
- Adjustable height settings
- Vinyl contoured hand grip for added comfort and long lasting wear

<b>PRODUCTS</b>	<b>Description</b>
WKSAN2B	STEEL ADULT WALKER 2 BUTTON 4EA/CV