

ProBasics Folding Walkcane Hemi Walker, 4/CA

HCPCS:

- The ProBasics Folding Walkcane Hemi Walker provides patients with the comfortable, stabilized travel assistance they need.
- Its aluminum frame allows it to be durable while still being lightweight and compact.
- Durable, anodized aluminum frame
- Lightweight compact design folds easily for easy storage and transport
- Special angled foot pieces provide stability
- Includes handgrips for added convenience and comfort
- Adjustable height from 28.5" - 33.5" to accommodate a range of patients
- 2-year warranty on frame
- Quantity: 4/case
- For comfortable and durable walking assistance, purchase the ProBasics Folding Walkcane Hemi Walker.

PRODUCTS	Description
PB6055	FOLDING WALKCANE HEMI WALKER 4EA/CV