

ProBasics Heavy Duty K0007 Wheelchair, 22" x 18" Seat with Footrests,

HCPCS:



- The ProBasics K7 Heavy Duty Wheelchair features a robust design to accommodate patients with wider seat and higher weight capacity requirements.
- With seat widths from 24 to 28 inches, a weight capacity of 450 - 600 pounds and a heavy-duty slide tube frame, the ProBasics K7 offers a rugged yet attractive option for the bariatric patient population.
- Heavy Duty, slide tube silver vein steel frame
- 8" x 2" semi-pneumatic casters are adjustable in 3 positions
- 24" rear urethane wheels with durable hand rims
- Removable desk length arms with padded armrests
- Heavy duty vinyl padded upholstery
- Easy push-to-lock wheel locks
- HCPCS* Code: K0007, K00195 or E0990 (elevating legrests only)
- WARNING: This product can expose you to chemicals including Lead, which is known to the State of California to cause cancer or birth defects or other reproductive harm.
- For more information, visit www.P65Warnings.ca.gov .

PRODUCTS

Description

WC72218DS

K7 WHEELCHAIR, 22X18 W/ SAFR 1EA/