

## ProBasics K2 Wheelchair with 16" x 16" Seat and Elevating Legrests

HCPCS:

- The ProBasics K2 Standard Hemi Wheelchair is designed to meet a variety of patient profiles.
- With a patient weight capacity of 300 pounds, a dual axle for hemi-height adjustment and multi-position caster forks, it is one of the more versatile wheelchairs on the market today.
- Standard features include flip back removable padded armrests, easy-to-clean vinyl upholstery, chart pocket and heel loops on the front riggings.
- Slide tube, powder-coated silver vein steel frame
- Dual axle provides easy transition of seat height to hemi-level
- Easy-to-clean, comfortable vinyl upholstery
- Flip-back, padded, upholstered armrests are removable
- Chart pocket on back for additional convenience
- Includes elevating padded legrests with composite footplates
- Elevating legrests are equipped with a pneumatic piston for infinite elevating positions
- Easy push-to-lock wheel locks

### PRODUCTS

### Description

WC21616DE

K2 WHEELCHAIR, 16X16 W/ ELR 1EA/