

ProBasics Offset Folding Cane, Black

HCPCS:

- ProBasics Folding Canes fold easily for storage or travel.
- The comfortable T-shape grip reduces hand stress and fatigue on the wrist for maximum comfort.
- Available in two attractive finishes with a limited lifetime warranty.
- Supports patient weights up to 300 pounds.
- Folds in four convenient parts for easy storage
- Comfortable, contoured T-shape handle
- Easy-to-use, one button height adjustment with a locking anti-rattle ring allowing for quick safe fitting; adjusts in 1" increments.
- Attractive powder-coated finish in black
- 300 lb weight capacity
- HCPCS* Code: E0100

PRODUCTS	Description
CNFBK	FOLDING CANE, BLACK 6/CS 6EA/CV