

ProBasics Sure Lever Release Folding Walker with 5" Wheels, (Adult)

HCPCS:

- The ProBasics Sure Lever Release Walker is an ideal design for patients with manual dexterity limitations.
- The sure lever release folding mechanism allows users to close the walker without removing their hands from the walker handle.
- Ample height adjustment options coupled with a sure lever release folding mechanism make it easy to use and accommodating to the most demanding patient needs.
- Supports patient weights up to 350 pounds.
- Attractive, lightweight aluminum frame
- Sure lever release folding mechanism allows users to close the walker without removing their hands from the walker handle
- Each side operates independently to allow easy movement through narrow spaces and greater stability while standing
- Front cross brace near top of walker allows for full stride
- Adjustable height settings
- Vinyl contoured hand grip for added comfort and long lasting wear
- Rear glide cap allows for walker to slide easily and smoothly over most surfaces (come standard on models with 5" wheels)

PRODUCTS

Description

WKAAWSL

ALUM ADULT WALKER-TRIG W/WHLS 4EA/CV