

Nitro Sprint Cane Holder

HCPCS:



FIGURE A



FIGURE B



FIGURE C



FIGURE D



FIGURE E

The Nitro Cane Holder is designed to fit on the Nitro Sprint Rollator with built-in accessory mount to provide convenience at every turn

- Fits on all Nitro Sprint Rollators with built-in accessory mount
- Securely holds most standard canes on left side of rollator
- Easily mount cane or remove cane from rollator for seamless transition from cane to rollator

PRODUCTS

Description

102662-CANE	Nitro Cane Holder, 1BX
-------------	------------------------

SPECIFICATIONS

Description

Product Dimensions	3.8" (L) x 2.4" (W) x 2" (H)
--------------------	------------------------------

CA Prop 65	no
------------	----