

Viverity Offset Cane Replacement Hand Grip, Black

HCPCS:

- This replacement hand grip provides extra cushion and comfort to the Viverity Aluminum Offset Cane.
- It comes equipped with a wrist strap and is available in a variety of fashionable and fun colors.
- Simple replacement for aluminum offset canes
- Cardio-grips add extra cushion and comfort
- Easy to install
- Black; available in other fashionable and fun colors
- Sold individually

PRODUCTS

Description

M206-6BK	Offset Cane Replacement Hand 6EA/PKG
----------	--------------------------------------