

Universal Walker Tray with Cup Holder

HCPCS:



FIGURE A




FIGURE B



- Allows personal items to be carried from room to room
- Contains one cup holder
- Fits most manufacturers' walkers
- Made of easy-to-clean durable plastic
- Easy-to-install
- Size: 23" (L) x 17" (W) x 1.5" (H)

PRODUCTS

Description

10124	Walker Tray, Univ, 3/cv
-------	-------------------------

SPECIFICATIONS

Description

Warranty	Limited 1 Year
Dimensions	23" (L) x 17" (W) x 1.5" (H)