

Walker Accessory Bag

HCPCS:



Take your items with you wherever you go using the functional and versatile bag that easily secures to most walkers. The front and rear pockets are perfect for keys or a cell phone, and the large lined interior pocket with O-ring zipper closure holds all your essential items. The bag is constructed of lightweight nylon material that will remain sturdy even with extensive use.

- Adjustable bag loop fasteners attach to wide range of mobility aids
- Spacious interior pocket with zipper closure and key ring
- Three front pockets hold essentials such as keys, snacks or cell phone
- Two back pockets allow items to be accessed from users side
- Lined interior and pockets
- Treated with water resistant technology
- Spot clean with wet cloth

PRODUCTS	Description
RTL10256NV	Univ MobilityBag, Med Flat, Navy, 1EA
RTL10256PK	Univ MobilityBag, Med Flat, Pink, 1EA
SPECIFICATIONS	Description
Weight Capacity	6 lb
Dimensions	16" (L) x 9" (W) x 2.5" (D)
Retail Box Dimensions	11" (L) x 6.1" (W) x 1.4" (D)
Weight	0.5 lb