

## Bariatric Aluminum Folding Walker- Two Button

HCPCS:



- For Bariatric use
- Steel legs and side braces
- Precision design provides additional strength while adding minimal weight
- Wider and deeper frame design accommodates individuals up to 600 lb
- Easy push-button mechanisms may be operated by fingers, palms or side of hand
- Each side operates independently to allow easy movement through narrow spaces and greater stability while standing
- Vinyl-contoured hand grip

| PRODUCTS | Description                |
|----------|----------------------------|
| 10220-1  | Walker,Oversized,Adt.,1/ea |

  

| SPECIFICATIONS         | Description      |
|------------------------|------------------|
| Warranty               | Limited Lifetime |
| Handle Height          | 32" - 39"        |
| Weight Capacity        | 600 lb           |
| Carton Shipping Weight | 12 lb            |
| Overall Width          | 23"              |
| Weight                 | 9 lb             |