

## AirgoOne Touch Button Folding Walker



HCPCS:



FIGURE A



FIGURE B



FIGURE C



FIGURE D



FIGURE E



FIGURE F



FIGURE G

- One-button folding mechanism makes opening and closing the walker easy
- Lightweight and easy to maneuver
- Front bar also acts as a support rail for use over a commode bowl
- Height adjustable for a handle height of 31.5 to 38
- Includes plastic glide tips

**PRODUCTS**

**Description**

770-165

Airgo F.Walker, Tips, Small, Silver, 1E

**SPECIFICATIONS**

**Description**

Weight Capacity

300 lb