

Folding Walker, Two Button with 5" Wheels

HCPCS: E0143



FIGURE A



- Easy push-button release mechanisms
- Each side operates independently to allow easy movement through narrow spaces and greater stability while standing
- Sturdy 1" diameter construction ensures maximum strength
- Newly designed rear glide cap allows for walker to slide easily and smoothly over most surfaces (Figure A)
- Vinyl-contoured hand grip
- Comes with wheels (Replacement wheels item # 10109)
- Limited Lifetime Warranty

PRODUCTS	Description
10244-4	Walker,Steel, Adult 2Button 5inWhls,4/CS
SPECIFICATIONS	Description
Warranty	Limited Lifetime
Weight Capacity	350 lb
Height	32.5"-41.5"
Carton Shipping Weight	36.79 lb
Overall Width	22.5"
Weight	8 lb
Wheel Size	5"