



ProBasics Aluminum Transport Wheelchair, 19-inch, Blue

HCPCS:



FIGURE A



FIGURE B



FIGURE C

- The ProBasics Aluminum Transport Chair is a compact, lightweight design with an easy-to-maintain finish.
- The flared handles and handy back-release hinge allow the back to fold down flat and compact for easy storage and transport.
- Comes complete with seat belt, tool-free swing-away composite footrests and heel loops.
- Fits patient weights up to 300 pounds.
- Lightweight aluminum construction (only 20 lbs)
- Flared handles and handy back-release hinge allow the back to fold down flat and compact for easy storage and transport
- Stylish aluminum frame with permanent rear axel for strength and durability
- Color: blue
- Padded full-length permanent armrests
- Padded nylon upholstery is comfortable, lightweight and easy to clean
- Seat belt for added safety

PRODUCTS

Description

TCA1916BL

ALUM TRANSPORT CHAIR-BLUE 19X16 W/ SAFR