

ProBasics K1 Wheelchair with 18 x 16 Seat,

HCPCS:



- The ProBasics K1 Manual Wheelchair combines style and durability at an affordable price.
- Standard features include flip back removable padded armrests, easy-to-clean vinyl upholstery, chart pocket and heel loops on the front riggings.
- Supporting patient weights up to 250 pounds, the ProBasics K1 Standard Wheelchair is a popular choice among providers.
- Available with swing-away footrests or elevating legrests.
- Slide tube steel frame in black powder coat finish
- Single axle positioning
- Easy-to-clean, comfortable vinyl upholstery
- Flip-back, padded, upholstered armrests are removable
- Chart pocket on back for additional convenience
- Includes tool-free swing-away footrests with heel loops or elevating padded legrests with composite footplates
- Elevating legrests are equipped with a pneumatic piston for infinite elevating positions

PRODUCTS

Description

WC11816DE	K1 WHEELCHAIR, 18X16 W/ ELR 1EA/
-----------	----------------------------------