

## ProBasics K2 Wheelchair with 18" x 16" Seat and Swing-Away Footrests

HCPCS:



- The ProBasics K2 Standard Hemi Wheelchair is designed to meet a variety of patient profiles.
- With a patient weight capacity of 300 pounds, a dual axle for hemi-height adjustment and multi-position caster forks, it is one of the more versatile wheelchairs on the market today.
- Standard features include flip back removable padded armrests, easy-to-clean vinyl upholstery, chart pocket and heel loops on the front riggings.
- Slide tube, powder-coated silver vein steel frame
- Dual axle provides easy transition of seat height to hemi-level
- Easy-to-clean, comfortable vinyl upholstery
- Flip-back, padded, upholstered armrests are removable
- Chart pocket on back for additional convenience
- Includes tool-free swing-away footrests with heel loops
- Easy push-to-lock wheel locks
- Precision sealed wheel bearings in front and rear ensure long lasting performance and reliability

### PRODUCTS

### Description

WC21816DS

K2 WHEELCHAIR, 18X16 W/ SAFR 1EA/