

ProBasics Reclining Wheelchair, 16" x 17", Removable Desk Arms & ELRs

HCPCS:



- The ProBasics Reclining Manual Wheelchair provides full back and headrest support as patient reclines from 90 to 160 degrees.
- The 13-inch headrest and padded, full-length arms are removable.
- Black, easy-to-clean, comfortable vinyl upholstery that accommodates a wide range of users..
- The frame has a chip-resistant silver vein finish.
- Mag-style wheels are lightweight and set back on frame to prevent patient from tipping when in the reclined position.
- Slide tube steel frame with an attractive silver vein finish that is chip resistant and easy to maintain.
- Rear anti-tippers and Easy push-to-lock wheel locks come standard.
- Padded Desk-length, detachable armrests provide added patient comfort.
- Easy-to-clean, comfortable vinyl upholstery that accommodates a wide range of users.
- Chart pocket on backrest provides additional convenience.
- Headrest extension comes standard with head immobilizer.

PRODUCTS

Description

WCR1616E	RECLINING WHEELCHAIR, 16X16 W/ ELR
----------	------------------------------------