

## Deluxe Sentra Full Reclining Wheelchair



HCPCS: K0001 OR K0002+E1226+E0990 OR K0195,K0001 OR  
K0002+E1226+E0990 OR K0195



FIGURE A



FIGURE B



FIGURE C



FIGURE D



FIGURE E



FIGURE F



FIGURE G

- Hydraulic reclining mechanism allows for infinite adjustments up to 180°
- Durable, heavy-gauge naugahyde upholstery
- Carbon Steel frame with triple-coated chrome for an attractive, chip-proof, maintainable finish
- Composite Mag-style wheels with chrome hand rims are lightweight and maintenance free
- 8" front casters
- Padded armrests provide added patient comfort
- Wheels set back on frame prevent tipping
- Precision sealed wheel bearings in front and rear ensure long-lasting performance and reliability
- Rear anti-tippers standard
- Comes standard with swing-away elevating legrests
- Aluminum foot plates are attractive, crack-proof and lightweight
- Front caster forks are adjustable in two positions
- Carry pocket standard
- Headrest extension with cushioned head immobilizer standard
- Comes with push-to-lock wheel locks

PRODUCTS	Description
STD22RBDDA	Sentra,22"Recl.DetDskArm,1/cs
SPECIFICATIONS	Description
Warranty	Limited Lifetime
Weight Capacity	450 lbs
Seat To Floor Height	21"
Carton Shipping Weight	72 lbs
Back Height	23"J33"
CA Prop 65	[WARNING: This product can expose you to chemicals that are known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to <a href="http://www.P65Warnings.ca.gov">www.P65Warnings.ca.gov</a> .]
ACCESSORIES	Description
STDS1074	6"Hemi Caster,SentraHD,1ea, EA
STDS860	Forwrd Stablzrs/Univ.Fit ,1pr/bx, BX
STDS802	Anti-Tipprs,w/Wheels,Adj,1pr/bx, BX
STDS840	RearWheel,TreadedUrethane, Viper GT 1/ea, EA
STDS801	Brake Handle Ext.,8", 1/ea, EA
WASR	Limb Support, Right,1/cs, CS
WASL	Limb Support,Left,1/cs, CS

## ACCESSORIES

## Description

STDS831	Heel Loops, 1 pr, BX
STDS1034	Wheelchair Cane/Crutch Holder,1/ea, EA
STDS806	Anti-Fold/Theft Device for w/c,1/ea, BX
STDS834	OverheadAnti-Theft Bar,Single,1/ea, EA
STDS855	Seat Belt,60" Auto-type,1/ea, EA
STDS851	Seat Belt, Velcro Close,1/ea, EA
STDS850	Seat Belt, Clasp Type,1/ea, EA
STDS6008-1	Carry Bag,W/C,Oxygen,1/ea, EA