

Blue Streak Wheelchair

HCPCS: K0001-1,K0001-1

- Blue Powder-coated frame, black cross brace
- Flip-back, padded desk-length arm
- Nylon upholstery is durable, lightweight, attractive and easy to clean
- Solid Polyurethane tires mounted on composite wheels provide durability and low maintenance
- Comes with push-to-lock wheel locks
- Comes with calf strap (-SF models only)
- Comes with swing-away footrests or elevating leg rests

| PRODUCTS | Description |
|--------------|--|
| BLS16FBD-ELR | BlueStreak16"D/L FlipBackArms,Elev,1/ea |
| BLS16FBD-SF | BlueStreak16"D/L FlipBackArms,SwgFt,1/ea |
| BLS18FBD-ELR | BlueStreak18"D/L FlipBackArms,Elev,1/cs |
| BLS18FBD-SF | BlueStreak18"D/L FlipBackArms,SwgFt,1/cs |
| BLS20FBD-ELR | BlueStreak 20"D/FlipBackArms,EivLeg,1/ea |
| BLS20FBD-SF | BlueStreak 20"D/FlipBackArms,SwgFt,1/ea |

| SPECIFICATIONS | Description |
|------------------------|---|
| Warranty | Limited 3 Year |
| Product Weight | 41.25 lb |
| Weight Capacity | 300 lb |
| Seat To Floor Height | 19.5" |
| Carton Shipping Weight | 47 lb |
| Back Height | 16" |
| CA Prop 65 | [WARNING: This product can expose you to chemicals that are known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov .] |

| ACCESSORIES | Description |
|-------------|---|
| STDS833N | Rear Anti Tipper w.wheels SSP&CS,1/pr, PR |
| STDS850 | Seat Belt, Clasp Type,1/ea, EA |
| STDS851 | Seat Belt, Velcro Close,1/ea, EA |
| STDS855 | Seat Belt,60" Auto-type,1/ea, EA |
| STDS856 | Seat Belt,60"Velcro-type,1/ea, EA |
| STDS801GT | WheelLockExtension,ViperGT,SSP2,1/ea, EA |
| STDS806 | Anti-Fold/Theft Device for w/c,1/ea, BX |
| WASL | Limb Support,Left,1/cs, CS |

ACCESSORIES**Description**

| | |
|----------|--|
| WASR | Limb Support, Right,1/cs, CS |
| STDS1034 | Wheelchair Cane/Crutch Holder,1/ea, EA |
| STDS831 | Heel Loops, 1pr, BX |