

Bariatric Steel Transport Chair

HCPCS: E1039



FIGURE A



FIGURE B

- Reinforced steel frame provides added support
- Supports individuals up to 450 lb
- Dual, reinforced steel cross braces
- 12" rear flat free wheels
- Heavy-duty, nylon reinforced upholstery with a back carry pouch
- Removable, reversible desk-length arms
- Padded armrests provide additional comfort
- Companion activated wheel locks
- Available in blue or red frames
- Swing-away footrests included

PRODUCTS

Description

BTR20-B	TransportChr.HD,20"Blue,1/cs
BTR20-R	TransportChr.HD,20"Red,1/cs
BTR22-B	TransportChr.HD,Blue,22",1/cs
BTR22-R	TransportChr.HD,Red,22",1/cs

SPECIFICATIONS

Description

Warranty	Limited Lifetime
Seat Dimensions	20" (W) x 18" (D) 22" (W) x 18" (D)
Weight Capacity	450 lbs
Seat To Floor Height	19.5"
Carton Shipping Weight	50 lbs 54 lbs
Overall Width	24.5" (open) 9" (closed) 26.5" (open) 11.5" (closed)
Weight	45 lbs 49 lbs
Back Height	16"
CA Prop 65	[WARNING: This product can expose you to chemicals that are known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov .]