

TranSport Aluminum Transport Chair

HCPCS: E1038



FIGURE A



FIGURE B

- Weighs 17.5 lbs without footrests
- Aluminum frame is lightweight and strong
- Frame magnets keep footrests secure against frame to prevent swinging open when the chair is carried
- Padded flip-back arms are removable
- Upholstery is comfortable, lightweight and easy-to-clean
- Two-way back release allows to push or pull lever
- The transport chair conveniently folds for travel and storage (Figure A)
- Back folds forward
- Comes with tool-free height adjustable swing-away footrests
- Aluminum rear wheel locks
- Weight Capacity: 250 lbs

PRODUCTS	Description
TS19	Tran Sport Chair,FlipBack-RemArms,1/cs
SPECIFICATIONS	Description
Warranty	Limited Lifetime
Seat Dimensions	19" (W) x 16" (D)
Weight Capacity	250 lbs
Seat To Floor Height	18"
Carton Shipping Weight	19.5 lbs
Overall Width	24.5" (W) (open) 10" (W) (closed)
Back Height	16"