

Travelite Transport Chair

HCPCS:



FIGURE A



FIGURE B



FIGURE C



FIGURE D



FIGURE E



FIGURE F

- The Travelite transport chair is lightweight and sturdy aluminum transport chair
- It is easily transportable and comes with a compact carry bag with wheels making it ideal for traveling
- Folds for storage and transportation
- It is ideal for use indoors and outdoors
- 18" square nylon seat comes with lift-up straps
- Comes with wheel locks for 8" solid tires
- Comes with padded, flip-up desk arms and height adjustable footrests

PRODUCTS	Description
TC005GY	Travelite Chair in a Bag,1/ea
SPECIFICATIONS	Description
Warranty	Limited Lifetime
Product Weight	20.5 lb
Product Dimensions	33" (L) w/ riggings x 22.5" (W) x 37.5" (H)
Seat Dimensions	18" (W) x 18" (D)
Weight Capacity	250 lb
Seat To Floor Height	19"
Carton Shipping Weight	23.2 lb
Folded Dimensions	11" (L) x 10" (W) x 31.5" (H)
Back Height	17.5"
CA Prop 65	[WARNING: This product can expose you to chemicals that are known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov .]