

Precision 50 PSI EasyAir Compressor with Flow Gauge

HCPCS:

- The Precision EasyAir Compressor is a powerful, high performance air compressor that delivers at least 13 liters per minute at 50 PSI.
- It is easy to use and lightweight, providing patients with the most convenient compressor option.
- Powerful 1/5 horsepower compressor delivers a minimum of 13 L/min at 50 PSI
- Only 20.4 lbs, making it easy to carry and transport
- Automatic resetting thermal limit switch protects the unit from overheating
- Eliminates the need for pressure/flow calculations with a built-in flow gauge
- Durable and easy to maintain
- Includes J bracket, flow gauge and non-marring feet
- For a safe and lightweight yet powerful compressor option, purchase the Precision EasyAir Compressor.

PRODUCTS

Description

PM15F

PRECISION 50 PSI AIR COMPRESS 1EA/
