

SOCCEBALL PEDIATRIC NEBULIZER

HCPCS:

- Designed to reduce anxiety during treatments and encourage compliance, the child-friendly Roscoe Soccer Compressor Nebulizer provides effective aerosol treatments with a state-of-the-art piston compressor.
- Delivers 29 PSI at maximum pressure and 8 to 13 PSI operating pressure (liter flow range 5 - 8 LPM).
- Equipped with a convenient nebulizer parking post.
- Comes complete with compressor, disposable nebulizer, mask, mouthpiece, 7-foot kink-resistant tubing, five extra filters and a five year warranty.
- Also available with a reusable neb kit option (Model No. SOCCER-TRU & SOCCER-TRUWB-PB).
- Child-friendly design helps to reduce anxiety during treatments
- Provides smooth, quiet treatments
- Equipped with a convenient nebulizer parking post
- Complete with compressor, nebulizer, child mask, angled mouthpiece, tubing, five extra filters and instruction manual
- TruNeb Reusable Neb Kit available for purchase with the compressor or alone

PRODUCTS

Description

SOCCE-TRUWB-PB

SOCCE BALL NEBULIZER W/ BAG&REUS 6EA/CV