

STRENGTHBRACE Universal Hinged Knee Brace

HCPCS:

- The STRENGTHBRACE Universal Hinged Knee Brace provides stability and support for chronic knee pain and instability.
- A unique wrap-around design fits a wide range of patient sizes.
- The lightweight, breathable neoprene material provides optimal compression while keeping patients comfortable and compliant.
- Medial and lateral stays provide excellent support and promote proper alignment while the hidden hinges can be set to provide the desired level of flexion/extension range of motion.
- Fits 12"- 23": measure around center of the knee with leg extended
- Enclosed tubular donut buttress
- Polycentric hinge provides flexion/extension and can be easily set for desired range of motion
- Hinge breaks at 0°, 10°, 20°, 40°, 60° and 90°
- Shortened knee brace design decreases the material bulk for increased comfort especially when worn for prolonged periods of time
- Can be worn on left or right knee PDAC assigned code: L1832 and L1833
- 1 year warranty

PRODUCTS

Description

BK1832-2	StrengthBrace Hinged Knee 25EA/CV
----------	-----------------------------------
