

ComboCare E-Stim and Ultrasound Combo Professional Device

HCPCS:

- NOTE: CLASS II - Rx Required.
- The ComboCare® is a professional two-channel device that combines clinical ultrasound therapy with five electrotherapy treatment waveforms: IFC, Premod, Russian, EMS and TENS waveform.
- Its advanced operating system ensures a complete range of therapy.
- The large easy-to-read LCD panel allows the clinician to quickly set treatment parameters.
- A dual frequency 5cm² applicator comes standard with the device.
- Its portable design lets you take it to your patient or keep it at the clinic.
- Its easy to use, and is up and running within minutes.
- Portable therapy - provide ultrasound and e-stim simultaneously, great for at the clinic or On the go
- Two channels of electrical stimulation output
- Large easy-to-read graphical LCD display with exact dial performance and programmable timer
- Ultrasound has 10% - 100% duty cycle selections

PRODUCTS

Description

DQ7844	ComboCare 3EA/CV
--------	------------------
