

## InTENSity CX4 Clinical Electrotherapy and Ultrasound System

HCPCS:

- The InTENSity CX4 is an advanced combination electrotherapy and ultrasound system offering the practitioner a wide range of treatment options in one user-friendly, ergonomic design at a fraction of what other systems would cost.
- Equipped with a color, touchscreen menu-driven interface, the InTENSity CX4 intuitively groups and displays clinical protocols, guiding clinicians step-by-step to ensure the ideal therapeutic treatment.
- Parameters can be easily selected and adjusted.
- Clinicians can choose to use preset clinical protocols or customize their own.
- The comprehensive library provides recommended electrode placement to ensure better treatment outcomes and improved clinic efficiency.
- Bilingual touchscreen user-friendly interface
- Dual frequency ultrasound at 1 MHz and 3 MHz
- Low BNR ratio of 5:1
- Multiple waveform electrotherapy: Interferential, Premodulated, TENS (Symmetrical Biphasic, Asymmetrical Biphasic, RAAS\*), NMS, Russian, High Volt, and Microcurrent
- Rapid Agonist Antagonist Sequencing (RAAS\*) replicates typical firing patterns of muscle contractions during functional movement for the best therapeutic outcome
- Independent intensity and parameter controls for each channel

### PRODUCTS

### Description

PRODUCTS	Description
DQ8001	THERATOUCH CX4 W/ THERAPY CART 1EA/