

Quattro 2.5 Professional Electrotherapy Device

HCPCS:

- The microprocessor controlled Quattro 2.5 is a professional four-channel electrotherapy device that provides IFC and Premodulated, Russian, EMS and TENS waveforms with two separate timers.
- Choose between several different amplitude modulation options.
- The interferential and pre-modulated modes offer frequency modulation as well as a static frequency option.
- The portable design lets you take it to your patient or keep it at the clinic.
- 4-channel electrotherapy with separate timers for simultaneous treatment for two treatment areas or two patients at once
- Large easy-to-read graphical LCD display with exact dial performance and two programmable timers
- Preset protocols as well as customized user defined protocols
- Waveforms include TENS, EMS, Premod (IF-2P), IFC (IF-4P) and Russian
- Constant current and constant voltage modes
- 10-presets per waveform totaling 50-presets in all (plus professional mode)
- Designed for tabletop, therapy cart or mobile use

PRODUCTS

Description

DQ8450

Quattro 2.5 3EA/CV