

## SoundCare Plus Professional Ultrasound Device

HCPCS:

- NOTE: CLASS II - Rx Required
- The SoundCare Plus is a professional ultrasound device, designed by clinicians.
- With this unique design, the SoundCare Plus comes equipped with two applicators, 1cm<sup>2</sup> and 5cm<sup>2</sup>, that are plugged into the device allowing the user to switch between applicators with a single button push.
- The SoundCare Plus has 20 presets and operates in 1 MHz and 3 MHz frequencies.
- The portable design lets you take it to your patient or keep it at the clinic.
- Its easy to use, and up and running within minutes.
- For the treatment of chronic and acute muscular pain
- Includes two applicators, 1cm<sup>2</sup> and 5cm<sup>2</sup>
- Fully functional 1 MHz and 3 MHz frequencies
- Simultaneously switch between applicators with the push of a single button
- 20 user-defined presets with pulsed and continuous therapy

### PRODUCTS

### Description

---

DQ9275	SoundCare Plus - Professional 3EA/CV
--------	--------------------------------------

---