

## TheraMED Professional Standard All-Terry Cover (10" x 28")

HCPCS:

- The TheraMED Professional Standard All-Terry Cover insulates heat to provide comfort and protect TheraMED Moist Heat Pack users from overexposure.
- It is designed with hook and loop closures to keep hot pack secure and prevent unwrapping or slipping.
- Constructed with soft terry cloth for optimum comfort and insulation; also available with foam layer (HPC1028-FF)
- Hook and loop closures keep hot pack secure and prevent unwrapping or slipping
- Machine washable
- Standard 10" x 28" size; available in a variety of sizes
- Ideal for both high-volume clinical use or at-home therapy
- One-year warranty

<b>PRODUCTS</b>	<b>Description</b>
HPC1028-TC	Standard All Terry Cover, 40EA/CV