

## TheraMED Professional Standard Foam-Filled Terry Cover (10" x 28")

HCPCS:

- The TheraMED Professional Standard All-Terry Cover insulates heat to provide comfort and protect TheraMED Moist Heat Pack users from overexposure.
- It is designed with hook and loop closures to keep hot pack secure and prevent unwrapping or slipping.
- Constructed with a polyurethane foam layer between two layers of comfortable terry for optimum insulation; also available with only terry (HPC1028-TC)
- Hook and loop closures keep hot pack secure and prevent unwrapping or slipping
- Machine washable
- Standard 10" x 28" size; available in a variety of sizes
- Ideal for both high-volume clinical use or at-home therapy
- One-year warranty

### PRODUCTS

### Description

---

HPC1028-FF	Standard Terry Cover, Foam 30EA/CV
------------	------------------------------------

---