

TheraTouch EX4 E-Stim Machine

HCPCS:

- The TheraTouch EX4 is an advanced four-channel electrotherapy system offering a wide range of treatment options in one user-friendly, ergonomic design.
- Equipped with a bilingual menu-driven touchscreen interface, the EX4 guides clinicians step-by-step through treatment protocols to ensure the ideal therapeutic treatment including electrode placement for better treatment outcomes.
- 4 channels of stimulation with the following waveforms:
 - Interferential
 - Premodulated
 - Biphasic
 - Asymmetrical
 - Symmetrical
 - RAAS
 - NMS
 - Russian

PRODUCTS

Description

DQ7001	THERATOUCH EX4 w/Therapy Cart 1EA/
--------	------------------------------------
