

## THERATOUCHE UX2

HCPCS:

- Equipped with a menu-driven, color touchscreen interface that intuitively groups and displays clinical protocols, guiding clinicians step-by-step to ensure the ideal therapeutic treatment.
- Appropriate parameter settings
- Ultrasound applicator recommendations
- Dual-holster design to accommodate both size applicators
- Ergonomic designed sound head provides for an excellent ERA
- A few simple key presses are all that is needed to quick-start a treatment.
- Parameters can be easily selected and adjusted.
- Clinicians can choose to use preset clinical protocols or customize their own.
- feature allows user defined protocols to be stored for future use.
- The UX2 offers a comprehensive library with treatment education for better treatment outcomes and improved clinic efficiency.

### PRODUCTS

### Description

---

|        |                      |
|--------|----------------------|
| DQ2000 | THERATOUCHE UX2 1EA/ |
|--------|----------------------|

---