

BESTCARE ELEC SIT-TO-STAND

HCPCS:

- SA182 SIT-TO-STAND BESTSTAND
- The Bestlift SA182 is an affordable care solution.
- This model is compact and easy to maneuver.
- Medicare reimbursable E0630*.
- PRODUCT FEATURES:
- 400lb capacity
- Full electric power
- Low-profile dual front casters (4.5")
- Adjustable knee pad
- Detachable foot plate

PRODUCTS

Description

SA182	BESTCARE DLX ELECTRIC STAND ASSIST LIFT
-------	---