

BESTCARE ELECTRIC LIFT

HCPCS:

- The BestLift PL182
- A perfect blend of form and function for patient transfers using 100% electric lift for those up to 400 lbs.
- With an overall base length of 44.9, maximum width of 42.9, and base height of 4.5, it is the perfect solution for safely transferring patients in small rooms with low beds.
- PRODUCT FEATURES:
 - 400lb capacity
 - Low-profile dual front castors (4.5")
 - New electronics with 24/7 capability
 - New hand operated base opening
 - Works with ProBasics slings

PRODUCTS

Description

PL182	BESTCARE DELUXE ELECTRIC PATIENT LIFT
-------	---------------------------------------
