

BESTCARE ELECTRIC PATIENT LIFT

HCPCS:

- The BestLift® PL228 Full Body Patient Lift (Drop Ship)
- A perfect blend of form and function for patient transfers.
- With an overall base length of 44.9, maximum width of 42.9, and base height of 4.5, it is the perfect solution for transferring bariatric patients in small rooms with low beds.
- The PL228 features:
 - 500lb capacity
 - Dual front castors
 - New electronics with 24/7 capability
 - New hand operated base opening

PRODUCTS	Description
PL228	BESTCARE PREMIUM ELECTRIC PATIENT LIFT