

Bestcare Electric Patient Lift

HCPCS:

- PL400E Patient Lift
- This battery-powered patient lift provides an economic option that is reliable and convenient for caregivers without sacrificing patient comfort.
- It has a compact design with a low base so caregivers can get as close to the patient as possible, making lifts incredibly easy and efficient.
- Easily adaptable for a wide range of challenges, with a six-point spreader bar
- Low 4.5" base easily slides under beds and chairs
- 3" front casters and 4" rear casters with locking brakes
- Battery powered operation system is low-cost yet effective
- Easy, tool-free assembly and disassembly
- Optional digital weigh scale and optional two-point spreader bar
- Weight capacity: 400 lbs

PRODUCTS	Description
PL400HE	BESTCARE GENESIS ELECTRIC PATIENT LIFT