

## BESTCARE HYDRAULIC PATIENT

HCPCS:

- The BestLift® PL400 Full Body Patient Lift
- One of the latest ProCare products.
- Ergonomically designed to reduce and eliminate back injuries among professional caregivers.
- Providing comfort and dignity for users.
- Main features:
  - 400lb capacity
  - Lift from floor and over high beds
  - 6-point spreader bar for extra comfort
  - Can be upgraded to a full electric lift with optional upgrade kit
  - Optional Scale

PRODUCTS	Description
PL400H	BESTCARE HYDRAULIC PATIENT LIFT