

## BESTCARE PL273 CONTROL BOX 1EA/

HCPCS:

PRODUCTS	Description
DSC700	BESTCARE DIGITAL PATIENT 1EA/
KIT-TA23-PL400	BESTCARE UPGRADE KIT-PL400 TO ELECTRIC
PB0036CLAMP	CLAMPS ONLY - PB0036/37 1PR/
PBPLSM001	PATIENT LIFT HYDRAULIC 1EA/
PBST	FLOOR MOUNT TRAPEZE 1EA/
PL182	BESTCARE DELUXE ELECTRIC PATIENT LIFT
PL228	BESTCARE PREMIUM ELECTRIC PATIENT LIFT
PL400EF	BESTCARE ELECTRIC FOLDING PATIENT LIFT
PL400H	BESTCARE HYDRAULIC PATIENT LIFT
PL400HE	BESTCARE GENESIS ELECTRIC PATIENT LIFT
PL500	BESTCARE PATIENT LIFT
PL600	BESTCARE ELEC BARIATRIC PATIENT LIFT
SA182	BESTCARE DLX ELECTRIC STAND ASSIST LIFT
SA182H	BESTCARE HYDRAULIC STAND ASSIST LIFT
SA400	BESTCARE ELECTRIC STAND ASSIST LIFT
SA400HE	BESTCARE ELECTRIC STAND ASSIST LIFT
SA600E	BESTCARE ELECTRIC BARI STAND ASSIST LIFT
SL-HC70001	BESTCARE HOYER COMPATIBLE SLING-LG
SL-HC70002	BESTCARE HOYER COMPATIBLE SLING-MED
SL-PM111	BESTCARE PEDIATRIC MESH SLING
SL-PP121	BESTCARE PEDIATRIC NYLON PADDED SLING
SL-SA669B	BESTCARE BARI STAND ASSIST SUPPORT STRAP
SL-TP332	BESTCARE SANI SLING-MED
SL-UHM832	BESTCARE UNIVERSAL MESH SLING W/HDRST-SM
SLUM-883	BESTCARE UNIVERSAM MESH SLING-LG
SL-UPH864	BESTCARE UNIVERSAL DLX SLING W/HDRST-XL
STA450	BESTCARE STANDING TRANSFER AID
WP-PL400EL-BAT	BESTCARE BATTERY FOR LEGACY ELEC LIFTS
WP-SA-PUMP	BESTCARE SA182H HYDRAUL. PUMP 1EA/
WP-TC12-CBOX	BESTCARE PL273 CONTROL BOX 1EA/