

Clear Plastic Wheelchair / Walker / Commode Bag, 1.5 mil

HCPCS:



- This Clear Plastic Wheelchair, Walker, and Commode Bag is ideal for homecare and long-term care environments.
- It meets the accreditation requirements by various accreditation organizations such as JCAHO, ACHC, and CHAPS by keeping contaminated equipment separate from clean, ready-to-deliver equipment.
- 30" x 20" x 35", available in a variety of additional sizes
- Clear plastic material
- Thickness: 1.5 mil
- Quantity: 200/roll
- Box size: 27" x 9" x 9"
- To best protect your patient's wheelchair, walker, or commode, purchase your Clear Plastic Equipment Bag.
- Compass Health Brands Offers a variety of equipment bags and covers for all your medical care needs.

PRODUCTS

Description

BAG-302035R

WCHAIR/WLKR/COMMODE BAG 1.5 MIL 30X20X35