

1 Drawer Bedside Cabinet



HCPCS:



FIGURE A



FIGURE B

- Tops feature High Pressure Laminate, Cabinet is Hi-Pressure Thermo fused Melamine. Drawers and doors feature solid wood edges
- Available in nine most popular wood grain finishes
- Drawer slides are epoxy coated and have 100 lbs rating (Figure A)
- European style hinges provide 110° accessible opening and self closing feature (Figure B)
- Ventilated back allows air flow throughout cabinet interior (Figure C)
- Warranty: Seven years

PRODUCTS

Description

MC-BDSC-11

Mill Creek 1 Door, 1 Drawer Bedside