

3 Drawer Chest

HCPCS:

- High Pressure wood grained Laminate Cabinet finish, Buried wood edges with High Pressure Laminate Top, Drawer front and door feature raised real wood detail, Solid wood feet _x000B_
- Available in six most popular wood grain finishes
- Drawer slides are epoxy coated and have 100 lbs rating (Figure A)
- Ventilated back allows air flow throughout cabinet interior (Figure B)
- Available in Arched and Flat Toe Kicks (Figure C)

| PRODUCTS | Description |
|-------------|----------------------------|
| RS-CHST-03 | Riverside Chest 3drwr |
| RS-CHST-03W | Riverside Chest 3drwr-wide |