

## 4 Drawer Chest

HCPCS:



- High Pressure wood grained Laminate Cabinet finish, Buried wood edges with High Pressure Laminate Top, Drawer front and door feature raised real wood detail, Solid wood feet
- Available in six most popular wood grain finishes
- Drawer slides are epoxy coated and have 100 lbs rating (Figure A)
- Ventilated back allows air flow throughout cabinet interior (Figure B)
- Available in Arched and Flat Toe Kicks (Figure C)
- Warranty: Seven years

| PRODUCTS   | Description           |
|------------|-----------------------|
| RS-CHST-04 | Riverside Chest 4drwr |