

## Oslo 1 Drawer Bedside Cabinet

HCPCS:



FIGURE A



FIGURE B

- Top High Pressure Thermo Fused Laminate
- Full Extension, Ball Bearing Drawer slides with rubber stops provide smooth operation, quiet closure
- Locking top drawer is a standard
- European style hinges provide 110° accessible opening
- Contemporary drawer handles with brushed nickel finish
- Ventilated back allows air flow throughout cabinet interior
- Available in three finish options (OSLO Cherry, OSLO Oak and OSLO Figured Mahogany)
- Specifications: 18" (D) x 21" (W) x 29.5" (H)

### PRODUCTS

### Description

NOBS1D1D-C	Oslo Bedside Cabinet, 1D/1D, Cherry
NOBS1D1D-FM	Oslo Bedside Cabinet, 1D/1D, FIG Mahogany
NOBS1D1D-O	Oslo Bedside Cabinet, 1D/1D, Oak