

Oslo 4 Drawer Chest

HCPCS:



- Top High Pressure Thermo Fused Laminate
- Full Extension, Ball Bearing Drawer slides with rubber stops provide smooth operation, quiet closure
- Locking top drawer is a standard
- Contemporary drawer handles with brushed nickel finish
- Ventilated back allows air flow throughout cabinet interior
- Available in three finish options (OSLO Cherry, OSLO Oak and OSLO Figured Mahogany)
- Specifications: 18" (D) x 37" (W) x 29.5" (H)

PRODUCTS

Description

NOC4D-C	Oslo 4 Drawer Chest, Cherry
NOC4D-FM	Oslo 4 Drawer Chest, FIG Mahogany
NOC4D-O	Oslo 4 Drawer Chest Oak