

3-Position Recliner, Bariatric

HCPCS:



- 3 positions: upright (Figure A), deep recline and elevated leg rest/footrest (Figure B)
- Retractable lock bar secures the chair in the desired position
- Self-adjusting suspension in the seat prevents bottoming out
- Side panels "pop off" for easy cleaning and maintenance
- Large, blow-molded tray locks in 5 positions and can be stored on the side of the chair
- 5" casters (2 with locks) make it easy to move and secure the recliner
- Wide, comfortable armrests are angled to provide support in any reclined position
- Neutral-colored, sturdy steel frame is easy to clean and maintain
- Moisture barrier on seat prevents seepage
- Comfortable built-in headrest
- Foam and vinyl meet California Technical Bulletin 117

PRODUCTS	Description
D574EW-CHAR	3 Position Extra Wide Charcoal
D574EW-CHOC	3 Position Extra Wide Chocolate
D574EW-TAN	3 Position Extra Wide Tan
SPECIFICATIONS	Description
Warranty	Limited Lifetime on Frame, 1 Year on Upholstery and Moving Parts
Seat Dimensions	20" (D) x 25" (W) x 21" (H)
Weight Capacity	500 lb
Carton Shipping Weight	127 lb
Width Between Arms	24.5"
CA Prop 65	[WARNING: This product can expose you to chemicals that are known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov .]