

## ProBasics Three-Position Recliner, Rose

HCPCS:



- The PROBASICS® three-position recliner offers both comfort and durability for day room or patient room use.
- Featuring full recline, elevated footrest with elevated and full upright positions, this recliner provides maximum flexibility for users.
- Need to move the recliner? The convenient push bar on the back and four rolling casters let you glide the chair where it's needed! For user safety, the two rear steering casters can be locked to secure the chair in place.
- FEATURES
- Hide-away footrest folds completely below chair to make rising to a standing position easier
- Vinyl seat and back padding to help reduce bottoming out
- Adjustable height settings
- Rust resistant rear locking caster forks
- Easy-clean, blow-molded side panels can be fast with no tools, just press the pin
- 5 position, locking molded trays easily slide ON and OFF with built-in side storage when not in use
- Convenient overlap tray with high lip completes the chair and the tabs let you move the tray freely

### PRODUCTS

### Description

PB6074RW

3 POSITION RECLINER - ROSE 1EA/