

## Homecare Mattress Extender 36" x 6" x 4"

HCPCS:

- The HDC-GAPCUSH\_6 is the quick solution to convert a standard 80" homecare mattress to 84" for taller users.
- Slip this extender between the mattress and the headboard.
- Combine this with the X1PBHBEXT that extends ProBasics full electric beds to 84".
- Dimensions: 36" wide x 6" high x 4" deep

<b>PRODUCTS</b>	<b>Description</b>
HDC-GAPCUSH_6	FOAM MATTRESS EXTENDER, 36X6X4 4EA/CV