

## HOSPICE COMPOSITE TOP OBТ

HCPCS:



- ProBasics Hospice Overbed Table with Composite Top
- Item #: PB3037C
- The ProBasics Hospice Overbed Table provides bedridden patients with a stable surface to continue independently carrying out activities that matter to them.
- The non-tilt includes a 30" composite top for easy cleaning and a long service life.
- Easily assembled with minimal hardware
- Swivel ball bearing casters for smooth movement
- 30" W x 15" D durable composite top provides plenty of space for daily activities
- The stable "H" base features a silver tone powder coat finish
- Manually locks at any set height to prevent the tray from accidental adjustment
- Set the height from 27" to 42" to accommodate most patients
- Supports up to 50 lbs.

**PRODUCTS****Description**

PB3037C

HOSPICE COMPOSITE TOP OBТ 1EA/