

Hydraulic Pump Only for Hydraulic Patient Lifts

HCPCS:

- This is the replacement hydraulic pump for the ProBasics Hydraulic Patient Lifts (PBPLSM001, PB0087 and PB0087X).
- It is ONLY the pump and does not include any other parts of the patient lift.
- Heavy duty round tube steel construction with a chip resistant hammertone finish
- Allows transfers of patients weighing up to 450 lbs
- Compatible only with ProBasics Hydraulic Patient Lifts (PBPLSM001, PB0087, PB0087X)
- One-year warranty
- Choose this Hydraulic Pump part when it needs replacing on your ProBasics Hydraulic Patient Lift.

PRODUCTS

Description

PB0086BVP

PBPLSM001/PB0087 HYDRAULIC 1EA/